Kontrollantigene Salmonella

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en

Intended use

The control antigens are used to check the agglutinability of the sifin Anti-Salmonella test reagents and for quality control when carrying out the slide agglutination test.

Principle of the test

Agglutination using an antigen-antigen-antibody reaction between specific *Salmonella* antibodies and control antigens of characterised *Salmonella* serovars.

Composition

The Salmonella control antigens are killed bacterial suspensions of characterised Salmonella serovars in phosphate-buffered saline (PBS) with 0.5% formalin.

REF	Available specificities	Antigen formula based on the White-Kauffmann-Le Minor scheme	
TS 1501	Control antigen Salmonella Paratyphi A-OH	2,12:a:[1,5]	Usually occurs as a monophasic variant.
TS 1502	Control antigen Salmonella Paratyphi B-OH	1,4,[5],12:b:1,2	
TS 1503	Control antigen Salmonella Paratyphi C-OH	6,7,[Vi]:c:1,5	Use TS 1507 for Vi
TS 1504	Control antigen Salmonella Typhi-OH	9,12, [Vi]:d-	Use TS 1507 for Vi
TS 1505	Control antigen Salmonella Typhimurium-OH	1,4,[5],12:i:1,2	Often occurs as a monophasic variant.
TS 1506	Control antigen Salmonella Enteritidis-OH	1,9,12:g,m:-	

TS 1507 Control antigen Salmonella Vi

The package size is 5 ml.

Form in which product is supplied, shelf-life and storage

The control antigens are provided in liquid, ready-to-use form. If stored at 2...8 °C both before and after opening, they may be used up to the date given on the label. After use, the bottles must be properly closed.

Warnings and precautions

As these materials contain formalin, contact with the skin and mucous membranes must be avoided! In case of contact, rinse with plenty of water.

Materials and equipment not supplied

Glass slides, stirring rods, physiological saline (NaCl 9,0 g/l), disposal containers for infectious material, sifin test reagents

Test material and methodology

The control antigens are used for quality control of the sifin Anti-Salmonella test reagents to determine the O, H and Vi antigens listed above using the slide agglutination test.

The control antigens must be carefully shaken prior to use.

One drop of each control antigen (about 25 µl) is mixed on a slide with one drop of the sifin test reagent (about 25 µl) which is required to determine the antigen formula of the control antigen.

The slide should be placed on a dark surface. The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control must be performed at the same time using physiological saline in place of the test reagent.

Evaluation of the agglutination result

Visible agglutination after the sample has been tilted back and forth less than 20 times. The negative control must still show a milky turbidity.

Note for the user: The control antigens are used in accordance with the above intended use exclusively for quality control of the serological typing with the slide agglutination test when using sifin test reagents.

Explanation of the symbols used

	Batch code (Lot)		Use by YYYY-MM (MM = end of month)
	Catalogue number		Temperature limitation
	In Vitro Diagnostic Medical Device		Consult instructions for use
	Slide agglutination		

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