

Ninhydrin 3.5%

REF CK9058

Intended Use

This test is used to determine the ability of bacterial strains to hydrolyze hippurate by the action of the enzyme hippurate hydrolase

Background

The hippurate hydrolysis test determines the ability of group B streptococci, as well as other bacteria to enzymatically hydrolyze sodium hippurate.

In addition to group B streptococci this method can also be used to identify *Gardnerella vaginalis* and *Campylobacter* spp.

Hippuric acid is hydrolysed by the enzyme hippuricase to glycine and benzoic acid. Ninhydrin evokes a 5 step reaction beginning with the deamination of glycine to form hydrindantin, carbon dioxide and ammonia. A condensation reaction occurs with the hydrindantin, ammonia and residual Ninhydrin to produce the final intense blue/black coloured complex.

Precautions

This product is for in-vitro diagnostic use and should be used by properly trained laboratory professionals. Universal precautions should be taken in the handling, processing and discarding of all materials used to perform the test. Do not use reagents after the expiration date shown on the product label has expired.

Methods

Use a fresh 18-24 hour culture as older cultures could be less metabolically active and results from these may be unreliable.

1. Prepare a dense suspension of the test strain (at least McFarland 4) in 0.25ml of saline
2. Add one hippurate hydrolysis Diatab to the suspension, close and incubate at 35-37C for a minimum of 4 hours or overnight.
3. Following incubation add 5 drops of Ninhydrin 3.5% to the tube, reincubate for no longer than 10 minutes and observe for a colour change.

For a positive reaction a dark blue/black coloured complex will form. With a negative reaction there should be no colour change, a faint colour change should be disregarded.

If results are inconclusive the test must be repeated using a longer incubation period for the hippurate hydrolysis, reincubation of an already completed test will not give a more conclusive result

Results

Positive Reaction – appearance of an intense blue or blue/black within 10 minutes.

Negative Reaction – light grey or no colour change.

Limitations

Plastic tubes should be used as glass may yield a positive result.

Incubation of the sample for more than 30 minutes following the addition of the Ninhydrin Reagent may yield false positive results.

Quality Control

A quality control should be undertaken daily or immediately prior to use.

Bacteria

Positive control-

Streptococcus agalactiae ATCC 12386

Negative control-

Streptococcus pyogenes ATCC 12344

Shelf Life & Storage

The expiry date, storage temperature (fridge) and storage conditions are indicated on the outer package label.

Materials provided

Each pack contains 5 x 3ml dropit bottles containing ready to use Ninhydrin 3.5% reagent.

Materials required but not provided

Sterile loops or needles

Sterile physiological saline

Small tubes

References

Barrow, G.I. & Feltham, R.K.A. Cowan and Steel's Manual for the Identification of Medical Bacteria. Third edition.

Ford, M. Medical Microbiology. Oxford University Press.

Jolly, J.L.S. Minimal criteria for the identification of *Gardnerella vaginalis* isolated from the vagina. J. Clin. Pathol. 1983;**36**:476-478.

Greenwood, J.R. & Pickett, M.J. Salient features of *Haemophilus vaginalis*. J. Clin. Microbio 1979;**9**:200-204.

Cacho, J.B. et al. Evaluation of a disk method for detection of hippurate hydrolysis by *Campylobacter* spp. 1989;**27**:359-360.

REF

Catalogue number

LOT

Batch number



Use by date

IVD

In-Vitro Diagnostic device



Contains sufficient for <n> tests



Temperature storage limitations



Consult instructions for use



Manufacturer

Issue	Date	Comments
3	29/10/2020	IFU format revision.

