

CE

40/41 Adeno-Strip

Rapid diagnostic test for *in vitro* detection of Adenovirus 40/41 in stool specimens

- **Easy to use**
- **Easy to read**
- **Specific**
- **Sensitive**
- **Reliable**
- **Rapid format (10')**




www.corisbio.com

CORIS BioConcept



Accurate, economical and fast procedure!



Performance compared to EIA (n = 153)

Specificity	**99.3% (143/144)**
Sensitivity	**100% (9/9)**
PPV	**90% (9/10)**
NPV	**100% (143/143)**
Accuracy	**99.4% (152/153)**

Bibliography

For more information:

www.corisbio.com/public/science

Or send an e-mail to: info@corisbio.com

Ordering Information

Enteric range (25 tests)

Description	Code
Rota-Strip	C-1001
Adeno-Strip	C-1002
40/41 Adeno-Strip	C-1003
Combi-Strip	C-1004
Crypto-Strip	C-1005
O157 Coli-Strip (20 tests)	C-1011
Giardia-Strip	C-1013
GastroVir-Strip	C-1016
Crypto/Giardia Duo-Strip	C-1018

Respiratory range (25 tests)

Description	Code
RSV Respi-Strip	C-1006
Adeno Respi-Strip	C-1009
Influ-A Respi-Strip	C-1010
Influ-A&B Respi-Strip	C-1012
Legionella V- TesT*	C-1815

All our products are also available in packaging of 10 tests, and in packaging of 20 tests, individually packed.
*Only available in 10 or 20 tests.

CORIS BioConcept

Science Park CREALYS
Rue Jean Sonet 4A - 5032 Gembloux - Belgium
Ph.: +32 (0)81 719 917 - Fx.: +32 (0)81 719 919
info@corisbio.com - www.corisbio.com

Distributed by: